

SAMPLE STAGE FOR HFRR

CANADA OFFICES:

Montreal: 514-937-7931 - 800-767-9695

Kitchener: 519-748-4612 - 800-683-6546

Calgary: 403-291-6685



www.folioinstruments.com



SAMPLE STAGE FOR HFRR

This sample stage allows the HFRR to be placed on a lab bench and raises the HFRR, so the weights are off the surface.

The feet are adjustable, so the unit can be leveled.

Specifications

- ✚ Stainless steel post
- ✚ 8 X 12 inches black anodized steel.

Price: CAN \$ 600.00